

Element RC



scil nr	Parameter
107093	Element RC Comprehensive Rotor
	Alb, Alp, Alt, Amy, BUN, Ca, Crea, GGT, Glu, Phos, TBil, TCho, TP
107094	Element RC Pre-Surgical Rotor
	Alb, Alp, Alt, BUN, CREA, Glu, TP
107095	Element RC Electrolyte Rotor
	Cal, Cl, K, Magn, Na, Phos, tCO2
107096	Element RC Equine/Large Animal Rotor
	Alb, Alp, Ast, BUN, Ca, CK, Crea, GGT, Glu, LDH, TBil, Trig, TP
108022	Element RC Kidney Rotor
	Alb, Amy, BUN, Ca, Crea, Glu, Phos, tCO2
108023	Element RC Bovine Monitoring Rotor
	Ast, Ca, CK, Magn, Phos
108024	Element RC Kidney Monitoring Rotor
	BUN, Crea, Phos
108025	Element RC Liver Rotor
	Alb, Alp, Alt, Ast, Chol, GGT, BA, TBil, TP
108085	Element RC Avian Reptilian Rotor
	Alb, Ast, Ca, CK, Glu, K, Na, Phos, BA, TP, UA
108668	Element RC Comprehensive Plus Rotor
	Alb, Alp, Alt, Amy, BUN, Ca, Crea, Glu, K, Na, Phos, TBil, TP