

# Fuji DRI-CHEM IMMUNO AU10V

## Reference Values



Parameter	Unit	Measurement Range	Dog*	Cat*	
<b>T4</b>	µg/dl	0.5 – 8.0	1.3 – 2.9	0.6 – 3.9	
	nmol/l	6.4 – 103.0	16.7 – 37.3	7.72 – 50.19	
<b>TSH</b>	ng/ml	0.25 - 5	< 0.5	-	
<b>Cortisol</b>	µg/dl	1.0 – 30.0 (Dilution: 1.0 - 50.0)	1.0 - 6.0 **	-	
	nmol/l	27.6 - 828.0 (Dilution: 27.6 - 1380.0)	28 - 166 **	-	
<b>Bile Acid</b>	µmol/l	2.0 – 150.0	fasting: < 7.9 postprandial: < 24.5	fasting: < 5 postprandial: < 15	
<b>Progesterone</b>	ng/ml	0.2 – 40	0.4 - 1.5		
	nmol/l	0.6 – 127.2	1.3 - 4.6		
			3.7 - 7.9 11.7 - 25.0		<b>Basal</b> (value between ovulations)
			11.7 - 28.2 37.3 - 89.8		<b>Ovulation</b> (value during ovulation)
			1.5 - 3.2 4.9 - 10.2	<b>Pregnancy</b> (value in pregnancy)	
				<b>Pre-Partum</b> (value before birth)	
<b>SAA</b>	µg/ml	3,75 – 225,0		0 – 5,49	
	mg/l	3,75 – 225,0		0 – 5,49	

### References:

\* Manufacturer´s Data

\*\* ACTH Stimulationstest. Low/High Dose Dexamethason Test



veterinary excellence

scil animal care company GmbH  
 Dina-Weissmann-Allee 6  
 68519 Viernheim  
 Tel.: +49 (0) 6204 78 90 222  
 Fax: +49 (0) 6204 78 90 200  
 Email: info-de@scilvet.com  
 www.scilvet.com