

# Fuji DRI-CHEM

## Reference Values

Parameter	SI Unit	Turtle	African Gray Parrot*	Bearded Dragon***	Iguana*	Boa**
ALB	g/l	13 - 30	15.7 - 32.3		13 - 16	25 - 39
ALP	U/l	36 - 156	20 - 160			
AMYL	U/l		210 - 530			
BIL	µmol/l	1.71 - 10.26			6.84 - 17.10	3.42 - 6.84
BUN	mmol/l	6.76 - 11.75	1.06 - 1.92	0.36 - 0.71	2.13 - 5.34	< 3.56
CA	mmol/l	2.49 - 3.61	2.11 - 3.24	2.14 - 6.78	2.24 - 6.25	2.24 - 5.48
CK	U/l		165 - 412			
Cl	mmol/l			80 - 140		
CREA	µmol/l	8.8 - 35.36	8.8 - 35.36		8.8 - 61.88	< 26.52
GGT	U/l	5 - 20	1 - 10			
GLU	mmol/l		10.54 - 19.42	7.71 - 16.15	8.32 - 15.54	0.55 - 3.88
GOT/AST	U/l	14 - 18	100 - 365	4 - 40		
GPT/ALT	U/l		5.0 - 12.0	3.0 - 5.0		
IP	mmol/l		1.03 - 1.74	1.13 - 3.16	1.13 - 3.16	0.83 - 1.58
K	mmol/l		2.9 - 4.6	1.0 - 6.5		
NA	mmol/l		157 - 165	141 - 190		
TCHO	mmol/l		4.14 - 11.0	8.08 - 31.70	2.84 - 8.83	1.37 - 3.60
TP	g/l	30 - 70			28 - 69	46 - 80
UA	µmol/l		276.68 - 565.10	95.17 - 678.13	89.23 - 356.91	71.38 - 345.01

### Bibliography:

- \* Tully, Dorrestein, Jones, Second Edition, Avian Medicin
- \*\* Avian Medicine and Surgery, Robert Altmann Reptile Medicine and Surgery, Douglas Mader
- \*\*\* Exotiv DVM Volume 8.5, Zoological Education Network, 2007

Vers. 20130827ENG



scil animal care company GmbH  
 Dina-Weissmann-Allee 6  
 68519 Viernheim  
 Tel.: +49 [0] 6204 78 90 222  
 Fax: +49 [0] 6204 78 90 200  
 Email: info-de@scilvet.com  
 www.scilvet.com