

QuickVet®

Reference Values



Parameter	Unit	Dog	Cat	Horse
PT	sec.	14 - 19	15 - 21	-
aPTT	sec.	75 - 105	94 - 125	-
equine Fibrinogen	g/l	-	-	1.5 - 4.0
canine Fibrinogen	g/l	1.2 - 3.0	-	-

References:

Manufacturer´s Data: Scandinavian Micro Biodevices ApS

Vers.: 20140623ENG



scil animal care company GmbH
 Dina-Weissmann-Allee 6
 68519 Viernheim
 Tel.: +49 (0) 6204 78 90 222
 Fax: +49 (0) 6204 78 90 200
 Email: info-de@scilvet.com
www.scilvet.com